OPTOSLIDE[®]

Disposable Counting Chamber

OPTOSLIDE are a disposable counting chamber to analyze motility and concentration of semen sample.



OPTOSCOPIA, Clever solutions for semen analysis



Optoslide

The OPTOSLIDE Counting Chamber's superior construction ensures a precise and consistent depth of the sample, facilitating accurate tracking of sperm movement and reliable estimation of concentration. With perfectly calibrated chambers, you can trust that each analysis will yield reproducible results.



OPTOSLIDE 204

Mat					■20 micron DPTOSLIDE * www.optoscopia.com
	1	2 3	4 5	6	0

OPTOSLIDE 206

Reference	OS20-4	OS20-6
Depth	20 micron	20 micron
Chambers	4	6
Sample	3μL	3µL
Shelf life	24 months	24 months
Application	Human/Animal/Industry	Human/Animal/Industry



Features



Our Optoslides are CE-marked fixed depth counting chambers used for both human and veterinary semen analysis. Manufactured in clean room, calibrated slides provide accuracy and reproducibility for semen analysis.

Optoslides are covered with special coating material to avoid semen for sticking to the slide and to avoid bubble formation. Sample filling and distribution by capillarity is guaranteed.





The material used for Optoslide is non-spermicidial, perfect for motility assessment.

Optoslides are single use counting chambers. This fact avoids time consuming cleaning procedures and minimize the risk of contamination related with reusable counting chambers.





Optoslide OS20-4 with 20 microns fixwed depth allow the spermatozoa to move freely for accurate motility assessment.

Optoslide Compatible with SIA lab and SIA lab plus and all CASA systems and microscopes.







MANUFACTURER

SPAIN :

calle Esmeralda 9, 08950 Esplugues de Llobregat Barcelona

+34 01 02 03 04 05

 TUNISIA:
43, rue de l'Energie,Zone Industrielle Charguia 1 Tunis.

L +216 01 02 03 04 05

🚖 contact@optoscopia.com

www.optoscopia.com

OPTOSCOPIA, Clever solutions for semen analysis

DISTRIBUTOR