Sperm Intelligent Analyser SIA lab®

Computer Assisted Semen Analysis (CASA)




As a leading manufacturer of CASA (Computer-Aided Semen Analysis) systems for semen analysis, we focus on combining optical technology and artificial intelligence (AI) to analyze sperm samples.



SIA lab®



measurement.



Mot Motility

Spermatozoa tracking is made at 50Hz frame rate to give precise trajectories.



Morphology

45MP camera gives all details of the spermatozoa for precise classification.



Vit Vitality

Determination of Live-Dead spermatozoa is accurate at more than 90%.




Recommended for samples having more than 1M/ml of round cell.



Frg DNA Fragmentation


The module of DNA Fragmentation is compatible with Sperm Chromatin Dispersion test.




WHO 6th edition

SIA lab and SIA lab plus follow recommendations of WHO 6th edition.






Perfectly calibrated, Optoslide® ensures reliable estimation of concentration and motility.



SIA lab® offers various customized and meaningful reports.

With SIA lab plus quality control module, reliable results are guaranteed.







Motility & Con	Automatic	Automatic	
Morphology	Semi automatic	Automatic	
Vitality	Semi automatic	Automatic	
DNA Frag	Semi automatic	Automatic	
Leukocyte	Semi automatic	Automatic	
Application	Human	Human	

